

SUSPEND EVERY INDIGO CARE CORDLESS DEVICE IN THE ROOM IN INCREMENTS OF 5 MINUTES

NO LEADS OR CABLES AROUND A RESIDENT

REDUCE THE RISK OF FALLS & POTENTIAL INJURIES

PORTABLE

EASY TO USE

PLUG & PLAY

INTEGRATES TO MOST NURSE CALL SYSTEMS VIA OUR RECEIVERS (RRU'S)

COST EFFECTIVE

KEY SPECIFICATIONS:

- *Australian Design*
- *C-Tick Approved*
- *Anti-Bacterial Silicon Buttons*
- *Raised, Easy to Press Buttons*
- *2 Buttons: White SUSPENSE, White CANCEL*
- *ID Encoded RF Transmitter*
- *Momentary Back Light*
- *Input Supply: Lithium Coin Battery CR 2032*
- *Screw Kit included*
- *Optional Double-Sided Tape*
- *Suits Standard Australian Electrical Boxes and Mounting Enclosures*

Suspense Point - Cordless

(Battery Operated)



TGA Registered ID 223902
RCM Approved
AS/NZS 61000.6.3.2012
Made in Australia

Part number: MED-RTU-SUSP

Indigo Care has designed a Cordless Suspense Point that will suspend all the Indigo Care's Cordless sensor mats and neck/wrist pendants temporarily so that the care-givers can monitor the residents and/or perform general housekeeping with ease, without setting off the nurse call system.

Needing only a Lithium Coin Battery (CR 2032) to power, it can be installed anywhere, but ideally next to the door entrance. Simply use the included screw kit or double-sided tape if preferred (not included).

Programming the Suspense Point to the Indigo Care's Cordless Receiver (RRU) is easy but must be done separately for the SUSPENSE and CANCEL button. See our Cordless Programming Instructions for complete details.

When the care-givers need to enter the room and suspend the Indigo Care's Cordless mats and pendants, simply press the SUSPENSE button. The button will light-up momentarily. This will suspend all Cordless mats and pendants in the room for 5 minutes (factory default setting).

You will notice that once suspended, the LED light on the Cordless Receiver will change from a green continuous light to a blinking amber light. This is an indication that the Cordless devices are in suspension mode.

Upon approaching the end of the 5 minutes suspension, the Cordless Receiver will give you a 10 second buzz warning. If you require more time, simply press the SUSPENSE button again, and the 5 minutes cycle starts over. If you leave the room before the 5 minutes lapses, press the CANCEL button upon your departure. Alternatively, it will automatically reset after 5 minutes. At the end of the suspended time, the Cordless Receiver will have a continuous green light, indicating that every Cordless mats and pendants device in the room has been reset and is again active.

Note: The Suspense Point does *NOT* de-activate other types of Corded devices (eg. Corded pendant) and cordless devices that are not designed by Indigo Care.

INDIGO CARE
(Horizon Electronics Company, T/as)
2/93 Rushdale Street, Knoxfield, VIC, 3180
PH +61 3 9763-2988
www.indigocare.com.au
ABN 36007302801

ORDER INFORMATION
Email your purchase orders to
sales@indigocare.com.au

