

# Deluxe Wet Sensor Monitor

REDUCE THE RISK OF  
INCONTINENCE  
ACCIDENTS

PORTABLE

EASY TO USE

EASY TO MOVE FROM  
ROOM TO ROOM

INTEGRATES TO MOST  
NURSE CALL SYSTEMS

PLUG & PLAY  
OPERATION

DURABLE & RELIABLE

COST EFFECTIVE



## PRODUCT CODE: MED-BM-01WS

(TGA Registration: 223902 - RCM Approved - AS/NZS61000.6.3.2012)

Indigo Care provides a Deluxe Wet Sensor Monitor that is designed to be used in conjunction with one of our Wet Sensor Mats or Pads. If the resident wets their bed or chair, the sensor strips in the mat will activate the monitor and immediately alert staff. We also supply a "Y" cable, as a separate item, for you to connect to your existing nurse call system. The monitor runs on 3 "AAA" batteries and has a low battery indicator light. It houses an inbuilt speaker for standalone operation, if required. There are three volume settings (Low, Medium, High) and three different tones from which to choose. A Delay Switch Setting allows three different timing trigger options, from 0 – 2 seconds. If needed, you can switch the sound off completely, by selecting "Mute".

### KEY SPECIFICATIONS:

- Part No. MED-BM-01WS
- Dimensions: 100mm X 70mm
- Connection: Curly Cord on RJ12
- Low Battery Detector
- Battery Operated 3 x AAA
- Optional: 6V DC Adaptor
- Warranty: 12 months

### PRODUCT ACCESSORIES:

- Wet Sensor Pad
- Wet Sensor Mat
- Y Cable

INDIGO CARE  
(Horizon Electronics Company, T/as)  
2/93 Rushdale Street, Knoxfield, VIC, 3180  
PH +61 3 9763-2988  
[www.indigocare.com.au](http://www.indigocare.com.au)  
ABN 36007302801

ORDER INFORMATION  
Email your purchase orders to  
[sales@indigocare.com.au](mailto:sales@indigocare.com.au)

