

REDUCE THE RISK OF
INCONTINENCE
ACCIDENTS

LOW PROFILE,
NON-SLIP MATERIAL

PORTABLE

EASY TO USE

EASY TO MOVE FROM
ROOM TO ROOM

WATERPROOF

WASHABLE

INTEGRATES TO MOST
NURSE CALL SYSTEMS

PLUG & PLAY OPERATION


DURABLE & RELIABLE

COST EFFECTIVE

LARGE WET SENSOR MAT
(381mm X 254mm)



Wet Sensor Mat



PRODUCT CODE: MED-WP-02

(TGA Registration: 223902 – RCM Approved – AS/NZS61000.6.3.2012)

Indigo Care provides a Wet Sensor Mat that is designed to immediately alert staff should a resident wet their bed or chair. This is most helpful to residents with incontinence issues. The mat works in conjunction with the Deluxe Wet Sensor Monitor. The Wet Sensor Mat is placed on the bed, directly under the bed linen. The Deluxe Wet Sensor Monitor is connected to the wet sensor mat and then, with a Y cable, the mat is connected to your existing nurse call system via the monitor. If the resident were to have an “accident”, the mat’s sensor strips will detect the moisture and the alarm will sound alerting the staff. The mat is waterproof and easy to wipe down and clean. A set of instructions are printed on the back of the pad.

KEY SPECIFICATIONS:

- **Dimensions: 381mm X 254mm**
- **Connection: Curly Cord on RJ11 (4P4C)**
- **Resistance <10k Ohm**
- **IP Rating : 66**
- **White & Green PVC, Thickness 0.20mm**
- **Cord Length 2 Meters**
- **Warranty: 12 months**

PRODUCT ACCESSORIES:

- **Y Cable**
- **Deluxe Wet Sensor Monitor (use 3 x AAA batteries-not included)**
- **6V DC Adaptor (optional)**

INDIGO CARE
(Horizon Electronics Company, T/as)
2/93 Rushdale Street, Knoxfield, VIC, 3180
PH +61 3 9763-2988
www.indigocare.com.au
ABN 36007302801

ORDER INFORMATION
Email your purchase orders to
sales@indigocare.com.au