

Mini Sensor Pad - Cordless



PRODUCT CODES:

RJ11 Mini Pad & Transmitter: MED-MS2X-WT

Full Kit: Receiver Unit, RJ11 Mini Pad & Transmitter: MED-MS2K-WT6(X), MED-MS2K-W3T(X)

Indigo Care provides the Mini Sensor Pad that can be used with our cordless range. The mini pad can be used as:

- A call pad for those with less dexterity
- An out of chair sensor
- An alarm for someone when they first sit on the toilet
- An alarm for someone getting up from the toilet

Dimensionally small (11 cm x 6 cm), this mini pad is great for anyone with limited movement. The pad can be placed around the wrist and pressure can be made by the palm of the other hand or it can be placed within reach of any part of the body that can apply enough pressure on the pad.

As a toilet sensor, it can notify when a resident sits down or gets up from the toilet. If required, it can also be an economical out of chair pad which will notify when a resident gets up from their chair.

KEY SPECIFICATIONS:

- Dimensions: 11cm x 6cm
- RJ11 Connection
- Velcro tabs to attach pad around wrist, toilet seat or bed frame
- White PVC thickness 0.06cm
- Cord length 1.42m
- Usage Record on back of pad
- Easy Clean
- Warranty: 12 months

Transmitter:

- Connection: RJ11 4P4C socket
- No External Antenna
- LED low battery indicator
- Transmits up to 30m (line of sight)
- Frequency: 433.92MHz
- Uses a 3V battery-CR2032
- Warranty: 12 months

CORDLESS RANGE:

- Cordless Receiver (RRU)
- Cordless Transmitter (RTU)
- Cordless Bed Sensor Mat
- Cordless Under Mattress Sensor Mat
- Cordless Floor Sensor Mat
- Cordless Cushion Sensor Mat
- Cordless Chair Mat
- Cordless Waterproof Neck Pendant
- Cordless Waterproof Wrist Pendant
- Cordless PIR Motion Sensor
- Cordless Suspense Point
- Cordless Doorway Sensor

ORDER INFORMATION

Email your purchase orders to
sales@indigocare.com.au

INDIGO CARE
(Horizon Electronics Company, T/as)
2/93 Rushdale Street, Knoxfield, VIC, 3180
PH +61 3 9763-2988
www.indigocare.com.au
ABN 36007302801

REDUCES THE RISK OF
FALLS & INJURIES

CAN BE USED AS A CALL
PAD FOR THOSE WITH
LESS DEXTERITY, AS AN
OUT OF CHAIR SENSOR,
OR AS AN ALARM
WHEN SITTING OR
RAISING FROM A
TOILET SEAT

USE SOLVENT WIPES TO
DISINFECTANT & CLEAN

PORTABLE

EASY TO USE

EASY TO MOVE FROM
ROOM TO ROOM

INTEGRATES TO MOST
NURSE CALL SYSTEMS

PLUG & PLAY
OPERATION

DURABLE & RELIABLE

COST EFFECTIVE

EXTREMELY SMALL PAD
(11cm X 6cm)

