

REDUCE THE RISK OF  
INCONTINENCE  
ACCIDENTS

LOW PROFILE,  
NON-SLIP MATERIAL

PORTABLE

EASY TO USE

EASY TO MOVE FROM  
ROOM TO ROOM

WATERPROOF

WASHABLE

INTEGRATES TO MOST  
NURSE CALL SYSTEMS

PLUG & PLAY  
OPERATION


DURABLE & RELIABLE

COST EFFECTIVE

LONG WET SENSOR MAT  
(760mm X 110mm)



# Wet Sensor Pad



**PRODUCT CODE: MED-WP-01**

**(TGA Registration: 223902 – RCM Approved – AS/NZS61000.6.3.2012)**

Indigo Care provides a Wet Sensor Pad that is designed to immediately alert staff should a resident wet their bed or chair. This is most helpful to residents with incontinence issues. The pad works in conjunction with the Deluxe Wet Sensor Monitor. The Wet Sensor Pad is placed on the bed, directly under the bed linen. The Deluxe Wet Sensor Monitor is connected to the wet sensor pad and then, with a Y cable, the pad is connected to your existing nurse call system via the monitor. If the resident were to have an “accident”, the pad’s sensor strips will detect the moisture and the alarm will sound alerting the staff. The pad is waterproof and easy to wipe down and clean. A set of instructions is printed on the back of the pad.

## KEY SPECIFICATIONS:

- **Dimensions: 760mm X 110mm**
- **Connection: Curly Cord on RJ11 (4P4C)**
- **Resistance <10k Ohm**
- **IP Rating : 66**
- **White & Green PVC, Thickness 0.20mm**
- **Cord Length 2 Meters**
- **Warranty: 12 months**

## PRODUCT ACCESSORIES:

- **Y Cable**
- **Deluxe Wet Sensor Monitor (requires 3 x AAA batteries, not included)**
- **AC Adaptor (optional)**

INDIGO CARE  
(Horizon Electronics Company, T/as)  
2/93 Rushdale Street, Knoxfield, VIC, 3180  
PH +61 3 9763-2988  
[www.indigocare.com.au](http://www.indigocare.com.au)  
ABN 36007302801

**ORDER INFORMATION**  
Email your purchase orders to  
[sales@indigocare.com.au](mailto:sales@indigocare.com.au)