### PORTABLE

## EASY TO USE

3 TYPES OF "Y" CABLE ADAPTORS, T1, T2, T3

USED WITH: FLOOR SENSOR MATS FLOOR CUSHION MATS BED SENSOR MATS BED SENSOR PADS UNDER MATTRESS MATS CHAIR SENSOR PADS WET SENSOR MATS WET SENSOR PADS

INTEGRATES TO MOST NURSE CALL SYSTEMS EXTRA 6.3MM SOCKET FOR NURSE CALL PENDANTS AND/OR BUZZERS

PLUG & PLAY OPERATION

**3 METERS IN LENGTH** 

MOLDED BOTH ENDS

**12 MONTH WARRANTY** 

# 'Y' Cable Adaptor – 6.3mm



# PRODUCT CODE: MED-Y-LEAD-T(X) (TGA Registration: 223902 – RCM Approved – AS/NZS61000.6.3.2012)

Indigo Care provides a special "Y" Cable Adaptor to be used with all our Fall Prevention sensor mats and pads. The "Y" cable comes as a standard item when you purchase a 6.3mm Floor Mat or 6.3mm Floor Cushion Mat and is also sold as a separate item when required. The "Y" Cable will allow you to connect the sensor mats and pads to your existing nurse call pendant system. The resident will then have a nurse call bell pendant and a mat/pad connected to the one nurse call point. There are three types of "Y" cables to choose from. This choice will depend on your existing nurse call system. The cord is 3 meters long which will allow for greater coverage as well.

#### TYPE 1: MED-Y-LEAD-T1

- TYPE 1 : Stereo Plug
- Suitable for 1 and 2 button nurse pendants with 6.3mm stereo plug tip common

#### **TYPE 2: MED-Y-LEAD-T2**

- TYPE 2: Mono Plug
- Suitable for 1 button nurse pendant with 6.3mm mono plug

#### TYPE 3: MED-Y-LEAD-T3

- TYPE 3: Stereo Plug
- Suitable for 1 and 2 button nurse pendants with 6.3mm stereo plug sleeve common

#### INDIGO CARE (Horizon Electronics Company, T/as) 2/93 Rushdale Street, Knoxfield, VIC, 3180 PH +61 3 9763-2988 www.indigocare.com.au ABN 36007302801

ORDER INFORMATION Email your purchase orders to sales@indigocare.com.au