

Pull Cord Pendant

MED-PCP-V2.0



■ Key Specifications:

- Call plate is ceiling mounted
- Mounting Plate – White & Plastic
- Reassurance Indicator – Red LED
- Operating Voltage – 12 to 24VDC
- Cord is 2.5 metres & Nylon
- Cord breaking strength – 10 kg
- Call Plate Dimensions:
 - 110 mm (L) x 67 mm (W) x 14 mm (D)
- Momentary Action Switch
- Protection against moisture & dust
- Optional: Variable Cord Length

Indigo Care has designed the Pull Cord Pendant with the elderly resident in mind. The call plate is mounted on the ceiling so the resident will have easy access to the pull cord to initiate a call. The cord length can be variable to suit different ceiling heights. There is a red LED for call activity. The switch is momentary action and is reset when the cord is released.

Part Number: MED-PCP-V2.0

- ***Easy to Use***
- ***Cost Effective***
- ***Integrates to most Nurse Call Systems***



INDIGO CARE
(Horizon Electronics Company, T/as)
2/93 Rushdale Street, Knoxfield, VIC, 3180
PH +61 3 9763-2988
www.horizone.com.au
ABN 36007302801

ORDER INFORMATION
Email your purchase orders to
sales@indigocare.com.au