

# 6.3mm Floor Cushion Sensor Mat



## PRODUCT CODES:

**6.3mm Floor Cushion Sensor Mat only: MED-FM-04**

**6.3mm Floor Cushion Sensor Mat & Y Cable (X = 1,2, or 3): MED-FM-04-6T(X)**

Indigo Care provides a New 6.3mm pressure sensitive Premium Floor Cushion “Crash” Sensor Mat that is designed for residents who are at risk of falling out of their bed, or those who require assistance in arising from the bed or a chair. The 6.3mm connection is durable and robust with a low profile ensuring the cable is secure and less vulnerable to breakages from external forces. Unlike other floor mats, our Premium Mats is a three-piece sensor mat made of 40mm foldable firm foam with a detachable 6.3mm “Y” cable to reduce trip hazard and damage to the electrical cord and provide an economical replacement cost if required. The 6.3mm Y-Cable is 3 meter long with a spare 6.3mm socket allowing the nurse call bell pendant to plug in as well so that either pressing the call button or standing on the mat will alert/trigger the nurse call system’s alarm.

Once the resident leaves the bed or chair and steps on the floor cushion, the sensor mat will activate the nurse call system and alert the staff. Our Floor Cushion Sensor Mats provide the ultimate in Fall Injury Protection as we have added an extra 40 mm of foam padding in case the resident falls while attempting to stand or walk. This mat is designed to absorb and cushion the resident if they fall over at full impact.

Our Premium Floor Sensor Mats are designed to work either with a cable (corded) or completely cordless with our Wireless Receiver and Transmitter.

## KEY SPECIFICATIONS:

- 6.3mm Connection
- Dimensions : 750mm X 1800mm X 40mm (29.5” X 70” X 1.5”)
- Silver Printing, Resistance < 150 Ohm
- Blue PVC Soft Foam Design
- 40 mm of Foam Padding
- 2-Fold Mat for Easy Storage
- Wired Connection: Detachable 6.3mm Y-Cable (included)
- Wireless Connection: RRU (Wireless Receiver) & RTU-FM (Reverse Polarity Wireless Transmitter). Purchased separately or as a Kit.
- Cord Length 3 Meters

## PRODUCT ACCESSORIES:

Included with the Floor Sensor Mat:

1 x 6.3mm - Y Cable: T1, T2, or T3

Wiring Configuration:

- T1-Tip & Ring
- T2-Mono
- T3-Tip & Sleeve

## ORDER INFORMATION

Email your purchase orders to  
[sales@indigocare.com.au](mailto:sales@indigocare.com.au)

INDIGO CARE  
(Horizon Electronics Company, T/as)  
2/93 Rusdale Street, Knoxfield, VIC, 3180  
PH +61 3 9763-2988  
[www.horizone.com.au](http://www.horizone.com.au)  
ABN 36007302801

6.3mm CABLE  
CONNECTION

PREMIUM OUT OF BED  
CUSHION MAT

THICK FIRM FOAM TO  
REDUCE, IF NOT  
ELIMINATE INJURIES

INCLUDE A 3-PIECE  
SENSOR MAT AND A  
DETACHABLE USB  
CABLE WITH  
ADDITIONAL SOCKET  
FOR PENDANT

FOLDABLE FOR EASY  
STORAGE

EASY TO MOVE FROM  
ROOM TO ROOM

WASHABLE

INTEGRATES TO MOST  
NURSE CALL SYSTEMS

PLUG & PLAY  
OPERATION

COST EFFECTIVE

LARGE PRESSURE MAT  
(750mm X 1800mm X 40mm)

